



Nasal non-invasive ventilation of neonates with nCPAP-PS

The breathing patterns of a very premature infant are frequently erratic and inadequate. Nasal CPAP effectively supports the breathing of preterm infants. Synchronized nasal intermittent positive pressure ventilation results in a decreased work of breathing, stabilizes the chestwall and reduces asynchronous motion between the chest and the abdomen and improves tidal volumes and minute ventilation.

Only a well coordinated product portfolio consisting of adapted mechanical ventilators, breathing circuits and a full range of nasal interface, keeps the promise to provide superior results for non invasive ventilation of neonates.

HAMILTON MEDICAL is offering the full range of accessories and consumables for successful nCPAP-PS therapy.

HAMILTON MEDICAL offers you all from one source.
More information: www.hamilton-medical.com


Accessories & Consumables

HAMILTON nCPAP-PS



HAMILTON nCPAP-PS generator

HAMILTON nCPAP-PS generator is a single use product for nCPAP therapy / non invasive ventilation therapy in treating neonates and premature infants in intensive care units together with a mechanical ventilator from HAMILTON MEDICAL.



Nasal Prongs

Single packed (10 pieces in each package) for single use in combination with nCPAP-PS generators. Available in 6 different sizes.



Nasal Masks

Single packed (10 pieces in each package) for single use in combination with HAMILTON nCPAP-PS generator. Available in three different sizes (small, medium and large). Anatomical form with 2x diaphragm.



Bonnets - single use

Single use bonnets for fixation of HAMILTON nCPAP-PS generator. One package includes 10 bonnets. 7 sizes available from XS to XXL. Material cotton and elasthan.



NEW
Feb
2011

XL Nasal Mask

Single packed (10 pieces in each package) for single use. Now available in XL size. Anatomical form with two diaphragms.

Neonatal nCPAP-PS consumables

REF 281975	HAMILTON nCPAP-PS Starter Kit (20 sets incl. mask, prongs and bonnets)
REF 281976	HAMILTON nCPAP-PS generator 20/box
REF 281977	Prong nCPAP-PS xsmall 10/packaging
REF 281978	Prong nCPAP-PS small 10/packaging
REF 281979	Prong nCPAP-PS medium 10/packaging
REF 281980	Prong nCPAP-PS large 10/packaging
REF 281981	Prong nCPAP-PS xlarge 10/packaging
REF 281982	Prong nCPAP-PS wide large 10/packaging
REF 281983	Mask nCPAP-PS small 10/packaging
REF 281984	Mask nCPAP-PS medium 10/packaging
REF 281985	Mask nCPAP-PS large 10/packaging
REF 281933 NEW	Mask nCPAP-PS xlarge 10/packaging
REF 281986	Bonnet nCPAP-PS xsmall 10/packaging
REF 281987	Bonnet nCPAP-PS small 10/packaging
REF 281988	Bonnet nCPAP-PS medium 10/packaging
REF 281989	Bonnet nCPAP-PS large 10/packaging
REF 281990	Bonnet nCPAP-PS xlarge 10/packaging
REF 281991	Bonnet nCPAP-PS xxlarge 10/packaging
REF 281992	Bonnet nCPAP-PS xxxlarge 10/packaging
REF 281994	nCPAP-PS measurement tape

Compatible with nCPAP-PS option for:

HAMILTON-C2
HAMILTON-G5
HAMILTON-S1



IntelliTrig optimized

Your contact: